

12 MONTH FACTORY GUARANTEE

Certificate of Warranty

Shins Bio Health Product (hereinafter to be called SBHP) warrants this heater to be free from defects for a period of twelve (12) months from the date of the original purchase subject to the following conditions. SBHP will replace (at its discretion) the faulty product free of charge (other than delivery costs) providing the product is used in accordance with the operating instructions and providing the customer returns the faulty product to SBHP at the customer's own expense. This warranty does not cover any damages due to accidents, punctures, plastic casing melting as a result of contact with hot metal (if a towel has not been used to protect the heater), microwave damage or any other misuse of the product by the customer. This warranty limits SBHP liability to the replacement of the product. Replacement parts are warranted for the portion of the original warranty period. Products returned for replacement after the warranty has expired will be subject to a fee for free parts and delivery. SBHP is not liable for damage of to or, any incidental or consequential damages caused by the use of this product. As you have taken the full advantage of the oral demonstration on the usage of this product, there is no cash refund. This warranty is valid only in Australia.

To validate the warranty or request for a mail order, detach and please return back to: **Shin Bio-Health Products Pty Ltd.**



Tel: (03) 9525 8573 Fax: (03) 9527 7985
 Mob: 0411 860 902
 Email: info@shinbio.com.au
 Website: www.shinbio.com.au
 Stall: Queen Victoria Market "K" Shed

SHIN BIO - SHIN CLINIC / BIO HEALTH PTY LTD.
 ACN: 094 375 863 | ABN: 23 094 375 863

HEAD OFFICE:
 628 Inkerman Road Caulfield North 3161 VIC
MANUFACTURING PLANT:
 39 Renver Road Clayton 3168 VIC Australia
POSTAL ADDRESS:
 PO Box 2392 Caulfield Junction 3161 VIC Australia

Warranty Card:

Customer Name:

Contact No.:

Address:

..... Postcode:

Purchased on...../...../.....

At



Tel: (03) 9525 8573 Fax: (03) 9527 7985
 Mob: 0411 860 902
 Email: info@shinbio.com.au
 Website: www.shinbio.com.au
 Stall: Queen Victoria Market "K" Shed

SHIN BIO - SHIN CLINIC / BIO HEALTH PTY LTD.
 ACN: 094 375 863 | ABN: 23 094 375 863

HEAD OFFICE:
 628 Inkerman Road Caulfield North 3161 VIC
MANUFACTURING PLANT:
 39 Renver Road Clayton 3168 VIC Australia
POSTAL ADDRESS:
 PO Box 2392 Caulfield Junction 3161 VIC Australia



HOW IT WORKS ??

The **Super Bio Health Heat Pack** contains water, sodium acetate and an aluminium disc. Clicking the aluminium disc forces a few liquid molecules to flip to a solid state. This triggers a chain reaction, causing the sodium acetate solution to supercool and crystallize – thereby producing heat.

placing it in a pot of boiling water (see Quick Instruction). We recommend you first place a plate or thick cloth on the bottom of the pot to protect the heat pack. Boil the heat pack until it appears clear of crystals. This should take 5-15 minutes. Using tongs, remove the heat pack from the water.

The Super Bio Health Heat Pack will generate heat for 30 minutes up to two hours, dependent on the size of the pack. For longer-lasting heat, massage the heat pack as this will disperse the heat and keep the heat pack soft and malleable. When the Super Bio Health Heat Pack becomes hard, you can reuse it by

Super Bio Health Heat Pack is portable which can give you instantaneous heat in any situation at a flex of button.

To use as a cold pack, store the Super Bio Health Heat Pack in the refrigerator for approximately 1 hour. Do not freeze the heat pack.

QUICK INSTRUCTION



1. Flex the metal disc inside the heat pack until crystals begin to form; the heat pack will harden as the liquid inside solidifies.



2. Massage the heat pack to evenly disperse the heat.



3. To reuse, place the heat pack in a pot of boiling water; place a plate on the bottom of the pot to protect the heat pack. Boil until the heat pack appears clear of crystals.



4. Use tongs to safely remove the heat pack from the pot. If the heat pack reactivates by itself, return to the boil for an extra 3-5 minutes.



5. The heat pack can be reused straight after boiling by wrapping the pack in a dry towel. When the heat pack cools, flex the metal disc to reactivate the heat pack.

MADE IN AUSTRALIA

**Always use tongs to safely immerse and remove the heat pack from boiling water.
 Made of Non-Toxic Materials.**

Heat Pack Range

The Super Bio-Health Heat Packs come in many convenient sizes.



Neck Heat Pack: 420mm x 130mm

The Neck Heat Pack comes in a versatile size that can be used to treat neck, shoulder, and back pain, abdominal cramps and period pain, aching knees and pulled hamstrings.



Large Heat Pack: 260mm x 180mm

The Large Heat Pack is recommended for use on the lower body; it can generate enough heat to treat lower back/lumbar pain, aching hips and period pain.



Sports Heat Pack: 260mm x 140mm

The Sports Heat Pack is ideal for relieving painful joints eg: wrists, elbows, knees and ankles. It is small enough to use as a hand or foot warmer and is excellent for treating tennis elbow and arthritis. Its small size is perfect for use by children and small adults.



Small Heat Pack:
125mm x 90mm

Commonly used as a hand warmer, the Small Heat Pack is designed to fit inside your pocket, desk, car glove-box or handbag. Its small portable size makes it perfect for use by children and small adults, and appeals to golfers, skiers, snowboarders, students, and anyone working outdoors.



Bio Ceramic Stone Belt-Back:
1350mm x 115mm - 190mm

Containing dehydrated ceramic stones, the Bio Ceramic Stone Belt comes in a lightweight velcro pouch that hugs the painful area - giving your body greater comfort. Used in conjunction with the Large Heat Pack, the Bio Ceramic Stone Belt conducts heat from the heat pack 40mm deeper into your skin.



Bio Ceramic Stone Belt-Sports:
550mm x 105mm - 150mm

Like the Bio Ceramic Stone Back Belt, the Bio Ceramic Stone Sports Belt maximises heat from the heat packs. Its smaller size is perfect for relieving painful joints eg: wrists, elbows, knees and ankles.

Cooling Eye Mask:

The Cooling Eye Mask can be used as a cold pack or a heat pack. To alleviate migraine, puffy eyes and headache, place the eye mask in the refrigerator for 30 minutes to an hour, then place over your eyes. This will halt inflammation of the sinus (sinusitis).



Warning: Heat Pack can be torn or penetrated. Please do not store with any sharp or underneath heavy objects. Not suitable for children age under 3 years old.

Applications

Sport & Leisure:

- :: Sporting stores
- :: Gun Shops
- :: Ski Resorts
- :: Scuba Dive Shops & Instructional Schools
- :: Boating Supply & Repair Shops
- :: High School Athletic Departments
- :: College Athletic Departments
- :: Sports Medicine

Therapeutic:

- :: Osteopaths :: Intravenous Therapy Clinics
- :: Veterinarians :: Ambulance Services
- :: Police and Fire Departments
- :: Emergency Medical Teams (EMT)
- :: Orthopaedic Doctors :: Home Health Services
- :: Dentists :: Massage Therapist :: Pharmacies
- :: Chiropractors :: Cosmetic Surgery
- :: Dermatologists :: Health and Diet Centres
- :: Naturopathic Physicians :: Occupational Safety
- :: Health Clinics :: Hospitals :: Physiotherapists

Mail Order:

Products	Size	Price	Qty
Neck Bio-Health Heat Pack	420mm x 130mm	A\$45.00	<input type="checkbox"/>
Large Bio-Health Heat Pack	260mm x 180mm	A\$35.00	<input type="checkbox"/>
Sports Bio-Health Heat Pack	260mm x 140mm	A\$30.00	<input type="checkbox"/>
Small Bio Health Heat Pack	125mm x 90mm	A\$12.00	<input type="checkbox"/>
Bio Ceramic Stone Belt - Back	1350mm x 115-190mm	A\$25.00	<input type="checkbox"/>
Bio Ceramic Stone Belt - Sports	550mm x 105-150mm	A\$20.00	<input type="checkbox"/>
Cooling Eye Mask	Average face size	A\$15.00	<input type="checkbox"/>
Pack A	Includes Neck, Large, Small Bio-Health Heat Packs and Bio Ceramic Stone Belt-Back	(Save \$17) A\$100.00	<input type="checkbox"/>
Pack B	Includes Neck, Large, Sports, Small Bio-Health Heat Packs, Cooling Eye Mask and Bio Ceramic Stone Belt-Back	(Save \$32) A\$140.00	<input type="checkbox"/>

*All prices include GST. Please add additional \$10.00 for delivery costs all around Australia.
**Only payments by Money Order or Credit Card will be accepted. No C.O.D's!

Customer Name: Contact No.:

Address: Postcode:

Type of credit card (Visa, Mastercard, Bankcard):

Credit Card No.: Exp:

Authorisation Signature:

Amount to be charged A\$

